Exp. 3

Aim: To synthesize Aspiriu

Apparatus Required: Measuring cylinders, beaker, conical flask, droppers, glass sod, filter perper, distilled water, and other general glasswere, Buchner funnel.

Chemicals Required: Salicyclic Acid, acetic anhydride, Conc. H2SOy, ethanol, mothanol, Felly,

Principle: Aspirin, also called Acetyl Sologolic Acid. the world, and has been termed as Wonder Drug of the century.

Its preparation is done by salicyclic acid and a cetic anhydride in the presence of a cid catalyst. i.e. H2SOy. In the seaction, phenolic hydroxyl geoup of salicyclic aced is esterified with a cetic onlyside to obtain Aspliin

After it has been prepared, impuraties may be there so filtration, then recrystallization and finally validation of purity by Fellz test is done.

Keactions!

Procedure:

(A) Synthesis;

1) Take 1.89 of salicyclic acid (Mol wto 138.12 glmol) into 150 ml conical flask.

and doctoynt as test and

- 2) Add 2-7 eq. of Acetic anhydride (Mol. wt. 102.08g/mol, Deneity 1.08g/ml) using a measuring cylinder to the salicyclic acid.
- 3.) Now add 5-6 deops of conc. H2504 and stir until all salicyclic acid is dissolved.
- 4) Allow the searction mixture to react for 15-20 nimules, lapte it un disturbed.
- After 15 mins, solid Aspiriu is formed.

 Add 50ml of distilled water & suish for 2 ming.

 and then filter using a Buchner funnel.
 - 6.) Solid Aspission can be collected from funnel, after filtration.

(B) Recrystallization;

- 1.) Dissolve the crude product in 7nil of ethanol and 15nil of distilled water in a beaker.
- 2.) put it on water bath to heat till you get a clear solution.
- 3.) Take it off from water both and keep it in the ice both for vecoy stallization. (for 10-15 mins)
- 4) Filtel again ming buchner funnel, and we get the pull white soled Aspiriu,
- 5) (alculate the percentage yield.
- 6.) Determine the melting point of acetylsalicy die acid.

(C) Validation of Pulity: (Fells Test)

This test is important to check presence /absence f phenol gop.

of phenol gop. since we have deve acylation of salicyclic acid
phenol has been converted to corresponding acetale

li) for salicyclic acid: it has phonol part so it will show the test.

1.) Take a pinch of salicyclic acid

2) dissolve in methernal

3) Add 1-2 deops of Fells, it will give a purple-violet coloniation, which shows presence of phenology.

is for Aspirin It should give - we test.

1) Take a swall amount of aspirir prepared
2) dissolve in methanol

3) Add 1-2 deops of Fells, colour will not change, It will remain yellow as the colour of Fells

Chemical Structures:

olulate the persent of gilla.

Purple Colour

Depending on the conc. of the reagent, more to trichelate vings maybe formed

Observation, & Calculations:

To calculate the amount of acetic anhydeide needed for 1.8 grams of salicy die acid.

→ We've been asked to add 2.7 g for 1.8 g salicyclic acid

Mal, wt. of salicyclic Acid = 138.129/mel

Mol. wt. of Acetic Anhydride = 102.08 gluwl

Density of - = 1.08 g/ml

(RAMoles of salicyclic acid = Moles of Acetic anhydride

2.7 x 1.8 = Moles (Acetic anhyderide)

Groams of Acetic anhydeide = 2.7 × 1.8 138:12 × 102.08

Volume of Acetic anhydride = 2.7 x 1.8 x 102.08 mL

= 3.33ml

-> To calculate Plesentinge Yield

Percentage yield = (Actual Yield / Theoretical Yield) X100

from chemical equation, we know that are male of sallcylic acid will give one male of a celyl salicyclic acid.

10, we started with 1 by salicyclic acid.

Moles of salicyclic acid = 1.8 = 0.013 moles

188:12

Hence we should get 0.013 moles of a cetyl salicydic acid

Mol. wt of acetyl solivitic acid= 180.158 gland

0.013 moles = 2.34g of a cetyl salicylic acid

Li) Theoretical yield

Actual yield = grams of product obtained by preforming reaction

Percentage Yield = 1.5 = 64.10%

Result: 1.) Asplitu was synthesized togsoon scalicyclic acid

- 2) Actual Yield = 1.5g
- 3) Percentage Yield = 64.10%
- 4) Melting pt, of Aspisiu = 134-136°C
- 5) Fellz test confirms the absence of phenolóc group in Aspirin prepared.

Precautions: 1. Dry Canical flask should be used

- 2. Measuring cylinder should be used for acetic anhydride
- 3. Conc. HSOy should be used with care
- 4. Care should be taken to isolate crystals of aspilin as fee as possible.