

Observations and Calculations

Standardization of Sodium Thiosulfate

S. No.	Volume of sodium thiosulphate used from the burette
1	9.4 mL
2	9.5 mL
3	9.5 mL

Given concentration of copper sulphate = 0.005 M

Equivalents of Copper = Equivalents of Thiosulphate

(Normality x Volume) of Copper = (Normality x Volume) of Thiosulphate

$0.005 \times 10 \text{ mL} = \text{Normality of Thiosulphate} \times \text{Volume of Thiosulphate (Burette Reading)}$

Suppose Burette Reading is 9.5 mL

Normality of Thiosulphate = $(0.005 \times 10)/9.5$

Normality of Thiosulphate = 0.0053 M

Observations and Calculations

Estimation of Iodine


S. No.	Volume of sodium thiosulphate used from the burette
1	2.4 mL
2	2.5 mL
3	2.5 mL



Consider the above reaction of iodate ions present in salt sample reacting with sodium thiosulphate

$$6 \times M_1 V_1 \text{ Iodate} = M_2 V_2 \text{ Thiosulphate}$$

$$6 \times M_1 \times 50 = 0.0053 \times \text{Burette Reading}$$

 Value obtained by standardizing sodium thiosulphate against a standard copper solution

Suppose Burette Reading is 2.5 mL

$$6 \times M_1 \times 50 = 0.0053 \times 2.5$$

$$\text{Molarity of Iodate solution} = 4.4 \times 10^{-6}$$

Observations and Calculations

Thus, the solution contains 4.4×10^{-6} moles of Iodine per Liter of solution

The solution contains $4.4 \times 10^{-6} \times 126.90$ grams Iodine per Liter of solution

Atomic Weight of Iodine = 126.90 g

Solution contains 5.6×10^{-3} grams Iodine per Liter of solution

= 5.6 mg of Iodine per Liter of solution

Concentration in ppm

To report in ppm

Suppose we started with 10 g of salt,

10 g salt was dissolved in 50 mL of water

In 50 ml of solution;

$(5.6 \text{ mg} \times 50) / 1000 = 0.28 \text{ mg}$ of Iodine

Thus, 0.28 mg of Iodine in 10 g of salt

ppm means parts per million; which is mg per Kg

$(0.28/10) \times 1000 = 28 \text{ ppm}$ of Iodine was present in our salt sample